

LOTUS Sharp Singularity™

"Water surrounds the lotus flower, but does not wet its petals"
-- Buddha --



The meniscus is anchored at the inner diameter of the LOTUS emitter because its surface is hydrophobic. This produces less evaporation, lower voltages, better ionization efficiency, and a more consistent spray.

FIT FOSSILIONTECH



More Peptide & Protein IDs, more consistently

LOTUS vs Regular 10 μm Sharp Singularity Emitters comparison*:


- 3.3% more peptide IDs,
- 5.5 fold reduction peptide ID variability

*Method: HeLa 50 ng injections, Ultimate 3000; Exploris 480, 75 μm -200mm C18 1.7 μm Waters Column. All emitters were conditioned by flushing two runs before analysis.


The LOTUS emitters lock the nano-electrospray meniscus inside the ID of the emitter. For the same emitter geometry, this results in a smaller meniscus, with a more stable anchorage.



LOTUS:



Regular:



The Sharp Singularity™
nESI Emitters

The Sharp Singularity™ Emitters

Improve the quality of your data
with robust and repeatable nano-electrospray ionization

Stable and repeatable signals require extremely tight tolerances at the microscopic scale.


FIT's unique micro-machining process produces:

- **Extremely sharp emitters**
- **Constant inner diameter (down to 10µm ID)**
- **Very tight tolerances**

Quality and traceability

We have full control of the production, QA/QC, and delivery process. All nano-ESI emitters are inspected and tested. To give you full control of your process, each emitter comes with a complete traceability and quality control report. With microscopy photos

Available geometries



Available emitter geometries:

Ref. 20-20	ID=20µm, OD=365µm, L=20cm, α= 7.5°
Ref. 20-10.5	ID=20µm, OD=365µm, L=10.5cm, α= 7.5°
Ref. 20-07	ID=20µm, OD=365µm, L=7cm, α= 7.5°
Ref. 20-06.25	ID=20µm, OD=365µm, L=6.25cm, α= 7.5°
Ref. 20-05	ID=20µm, OD=365µm, L=5cm, α= 7.5°
Ref. 20-02	ID=20µm, OD=365µm, L=2cm, α= 7.5°
Ref. 10-20	ID=10µm, OD=365µm, L=20cm, α= 7.5°
Ref. 10-10.5	ID=10µm, OD=365µm, L=10.5cm, α= 7.5°
Ref. 10-07	ID=10µm, OD=365µm, L=7cm, α= 7.5°
Ref. 10-06.25	ID=10µm, OD=365µm, L=6.25cm, α= 7.5°
Ref. 10-05	ID=10µm, OD=365µm, L=5cm, α= 7.5°
Ref. 10-02	ID=10µm, OD=365µm, L=2cm, α= 7.5°

Other geometries upon request:

Available IDs: 100µm, 75µm, 50µm, 40µm, 30µm, 20µm, 15µm, 10µm
Available lengths: from 2cm to 50 cm

Ion source /column– emitter compatibility table:

Ion Source model	Emitter recommended by Ion Source provider	Our references
Nanospray Flex™ - Thermo/Proxeon	≤ 20µm ID, 360µm OD, 10.5 cm long ¹	Ref. 20-10.5 Ref. 10-10.5
NanoFlow Spray™, NanoFlow Z-Spray™, Micromass NanoFlow™, NanoLock Spray™, - Waters	20µm ID, 360µm OD, 2.5" or 6.25 cm long ²	Ref. 20-06.25 Ref. 10-06.25
CaptiveSpray™ - Bruker	20µm ID, 360µm OD, 2 cm long ³	Ref. 20-02 Ref. 10-02
NanoSpray® III Ion Source - SCIEX	10-20µm ID, 360µm OD, 7 cm long ⁴	Ref. 20-07 Ref. 20-07
SUPER SESI™ - FIT	20µm ID, 360µm OD, 20 cm long	Ref. 20-20
Flex Ion connect for µPAC columns - Pharmalfluidics	20µm ID, 360µm OD, 5 cm long	Ref. 20-05 Ref. 10-05

1. Nanospray Flex Series Ion Source. User Guide - 60053-97127/ Revision B - Thermo Fisher Scientific.
2. Universal NanoFlow Sprayer. Installation and Maintenance Guide - 71500110107/ Revision C - Waters Corporation.
3. The impact II, a very high-resolution Quadrupole time-of-flight instrument (QTOF) for deep shotgun proteomics. Mol Cell Proteomics. Beck S. et al. 2015;14(7):2014–2029. doi: 10.1074/mcp.M114.047407.
4. NanoSpray® III Ion Source Consumables - D5037833 A - AB SCIEX.

Find pricing and more: <https://www.fossiliontech.com/nanoesi-emitters>